

# BLOCH EQUATION COMPUTATIONAL SIMULATION AND LOW FIELD STRENGTH NMR

Harvard College CHEM 10

Abel Bellows

# Content

1 Nuclear Spins

2 Bloch Equations

3 NMR: conceptual  
model

4 Computational  
Simulation

5 Physical  
Demonstration


6 Successes and  
Failures

# Basics of Spin

- Fundamental particles have "Spin."
- Spin can be measured in x,y, and z directions.
- The nuclei of atoms also have spin.

# Bloch Sphere

- Spin state can be represented as a linear combination of x, y, and z spin wavefunctions.
- Represented as vectors with head on a unit sphere.
- Bloch Sphere



Bloch Sphere

$$\frac{N_m}{N} = \exp\left(\frac{-E_m}{k_B T}\right) / \sum_{m=-I}^I \exp\left(\frac{-E_m}{k_B T}\right)$$

$$\approx \frac{1}{2I+1} \left(1 + \frac{m\hbar\gamma B_0}{k_B T}\right),$$

$$M_0 = \gamma\hbar \sum_{m=-I}^I m N_m$$

- The relative population of protons in spin state  $m$  ( $N_m$ ) is given by the Boltzmann distribution.
- At  $T = 398K$  spins are evenly distributed
  - *thermal  $E \gg \hbar m\gamma$*
- The net magnetization ( $M_0$ ) of the sample is the sum of the magnetic moments of the protons in the population
- $M_0$  is a vector on a Bloch Sphere.

$$\frac{d\mathbf{M}(t)}{dt} = \mathbf{M}(t) \times \gamma \mathbf{B}(t).$$

$$\omega_0 = -\gamma B_0.$$

- $\frac{d\mathbf{M}(t)}{dt}$  in external magnetic field  $\mathbf{B}$  is given by the torque applied to  $\mathbf{M}(t)$  by  $\mathbf{B}$ .
- If  $B(t) = B_0$   $\mathbf{M}$  precesses with frequency equal to  $\omega_0$

In a completely stable system,  $\mathbf{M}$  would precess around  $\mathbf{B}$  without decay.

Since "thermal equilibrium must be re-established," the Bloch Equations give us a way to model how  $\mathbf{M}(t)$  returns to equilibrium.

$$\frac{dM_x(t)}{dt} = \gamma \left[ \mathbf{M}(t) \times \mathbf{B}(t) \right]_x - R_2 M_x(t)$$

$$= \gamma \left[ M_y(t) B_z(t) - M_z(t) B_y(t) \right] - R_2 M_x(t),$$

$$\frac{dM_y(t)}{dt} = \gamma \left[ \mathbf{M}(t) \times \mathbf{B}(t) \right]_y - R_2 M_y(t)$$

$$= \gamma \left[ M_z(t) B_x(t) - M_x(t) B_z(t) \right] - R_2 M_y(t),$$

$$\frac{dM_z(t)}{dt} = \gamma \left[ \mathbf{M}(t) \times \mathbf{B}(t) \right]_z - R_1 \left[ M_z(t) - M_0 \right]$$


$$= \gamma \left[ M_x(t) B_y(t) - M_y(t) B_x(t) \right] - R_1 \left[ M_z(t) - M_0 \right],$$

$R_1$ : Spin-lattice relaxation rate constant (**how fast does the spin distribution return to Boltzmann**)

$R_2$ : Spin-spin relaxation rate constant (**How fast do the spins relax in the plane transverse to  $B$** )



$t = 3.0359999999997767 \text{ s}$



Bloch Procession


Green:  $M$   
Red:  $B_0$

# How can we use Bloch Precession to study molecules?




- $R1$  and  $R2$  are determined by the chemical environment around a nucleus. If we could measure the precession of  $M$ , we could learn about a nucleus' position in a molecule.

# NMR Machines

- A sample is placed in a strong B field of around 10T
- The effective B field is temporarily changed by applying an RF pulse
- The RF pulse knocks down  $M$
- $M$  processes around  $B$  according to the Bloch Equations.
- The changing emf is converted into an electric signal



# Fast Fourier Transform (FFT) of Free Induction Decay (FID) signal is the NMR spectrum. Its peaks reveal the chemical structure of a given molecule.



# Computational Bloch Equation/NMR Simulator

- Demonstrates behavior of processing  $M$  in NMR machine
- Configure all aspects of the device including  $B_{DEV}$ ,  $M$ ,  $T1$ ,  $T2$
- Send RF pulses with specified  $B_{\text{effective}}$
- Source code and detailed instructions at [https://github.com/abellows50/CHEM10\\_NMR\\_PROJECT](https://github.com/abellows50/CHEM10_NMR_PROJECT)

 README

## NMR Machine Simulator

### Purpose

To demonstrate the macroscopic behavior of  $M$  in an NMR machine and to make clear the workflow of NMR analysis, I constructed a simulator of Bloch Equations. The user can configure all aspects of the device including  $B_{DEV}$ ,  $M$ ,  $T1$ ,  $T2$ , and can send RF pulses with a specific applied magnetic field.

### Usage

#### Python Environment Configuration

To setup the python environment, first create your venv by running `python3 -m venv venv` then activate your venv by running `source venv/bin/activate`. Install all of the required packages by running `pip3 install -r requirements.txt`.

#### Start the program

To run the simulator run `python3 NMR_SIM.py`. You may want to use a screen split function so that you can see all 4 windows at one time.

#### Configure simulation

To configure environmental variables of the simulator, use the `config` command in the command console. You can configure the magnetic field of the device ( $B_{DEV}$ ),  $T1$  and  $T2$  ( $T1$  and  $T2$ ), and net starting magnetization ( $M$ ) with the following commands

1. `config B_DEV <Bx> <By> <Bz>`
2. `config M <Bx> <By> <Bz>`
3. `config T1 <value>`
4. `config T2 <value>`

#### Start the simulation

You can start the simulation by running `start`. You can pause the simulation by running `pause`.

#### Initiating radio pulse

You can initiate a radio pulse with the following commands: `pulse <duration>` or `pulse <duration> <Bx> <By> <Bz>`

#### Capture a photo of the simulation

You can capture a photo of the current frame by running `capture <filename>.png`

#### Fourier Transform to see the final NMR spectrum

Viewing the FID graph run `fft <t1> <t2>` to initiate an Fourier transform with the data between  $t=t1$  and  $t=t2$ .

 Harvard College CHEM 10

Bellows

# Flow Chart of Code




Figure 2





Figure 1



CHEM10\_NMR\_PROJECT — python3 NMR\_SIM.py — Python NMR\_SIM.py — 46x11

..llows50/pset6

-zsh ...

python3

```
starting at t=0.3020000000000002
lasting for t=0.05
```


```
NMR cmd >> pause
rf pulse shut off at t=0.35200000000000026
paused simulation...
```

```
NMR cmd >> fft 0.4 1.2
```

```
NMR cmd >> █
```


VisPy canvas

t = 1.3769999999999951 s





# Physical Attempt

Using Andy Nicol's Hackaday "Nuclear Magnetic Resonance for Everyone," I built a basic low field NMR machine.









Clockwise from  
top left: control  
Box, coil with  
water sample,  
amplifier and  
oscilloscope detector



Results with a magnetic stirrer under the coil to create a string rotating magnetic field.



Voltage vs Time graphs for low, medium, and high speeds on the magnetic stirrer.  
Clockwise from top left.




Water sample after a 12 second polarization at 4 Amps and 17 V.  
There is too much noise to see the FID signal or get a clear fft.

### Sources of noise:

- RF noise inducing current in coil
- Weak connections
- Crossed wires inducing current
- Static electricity buildup
- Random thermal fluctuations
- Faulty equipment

### Mitigation

- Increase magnetic field strength to increase signal strength
- Use faraday cage around apparatus
- Ground more effectively



## Works Cited

- JOHN CAVANAGH, WAYNE J. FAIRBROTHER, ARTHUR G. PALMER III, MARK RANCE, & NICHOLAS J. SKELTON. (n.d.). CLASSICAL NMR SPECTROSCOPY. Protein NMR Spectroscopy (Second Edition) Principles and Practice.
- Kaseman, D. (2013, October 2). Relaxation. Chemistry LibreTexts.  
[https://chem.libretexts.org/Bookshelves/Physical\\_and\\_Theoretical\\_Chemistry\\_Textbook\\_Maps/Supplemental\\_Modules\\_\(Physical\\_and\\_Theoretical\\_Chemistry\)/Spectroscopy/Magnetic\\_Resonance\\_Spectroscopies/Nuclear\\_Magnetic\\_Resonance/NMR\\_-\\_Theory/Relaxation](https://chem.libretexts.org/Bookshelves/Physical_and_Theoretical_Chemistry_Textbook_Maps/Supplemental_Modules_(Physical_and_Theoretical_Chemistry)/Spectroscopy/Magnetic_Resonance_Spectroscopies/Nuclear_Magnetic_Resonance/NMR_-_Theory/Relaxation)
- Liao, Dr. S., Dourmashkin, Dr. P., & Belcher, P. J. W. (n.d.). Faraday's Law of Induction. MIT Physics 8.02 Electricity and Magnetism.  
<https://web.mit.edu/8.02t/www/802TEAL3D/visualizations/coursenotes/index.htm>
- Nave. (n.d.-a). Nuclear Magnetic Resonance. Retrieved November 16, 2025, from <http://hyperphysics.phy-astr.gsu.edu/hbase/Nuclear/nmr.html#c1>
- Nave. (n.d.-b). Nuclear Spin Polarization. Retrieved November 16, 2025, from <http://hyperphysics.phy-astr.gsu.edu/hbase/Nuclear/spinpol.html#c1>
- Nave. (n.d.-c). Nuclear Spin Relaxation. Retrieved November 16, 2025, from <http://hyperphysics.phy-astr.gsu.edu/hbase/Nuclear/spinrel.html#c1>
- Newman, M. (n.d.). Bloch sphere | Prefetch. Retrieved December 8, 2025, from <https://prefetch.eu/know/concept/bloch-sphere/>
- NMR Basics for the absolute novice. (n.d.). JOEL Solutions for Innovation. Retrieved December 8, 2025, from <https://www.jeolusa.com/RESOURCES/Analytical-Instruments/NMR-Basics>
- WileyPLUS - NMR spectroscopy and nuclear spin. (n.d.). Wiley Plus. Retrieved December 8, 2025, from <https://www.johnwiley.net.au/highered/chemistry2e/wileyplus/basics/nmr/>

## Acknowledgements

- Khaled: Thank you for spending so many hours advising me in this project and for providing wisdom and support throughout the journey
- Andrew: Thank you for helping me clarify my ideas and realize what was possible for this project.
- Professor Park and Professor Xu: Thank you for showing how magical and beautiful Quantum Mechanics is and for inspiring all of this work.